
# SAE Dental-Produktkatalog

SAE Catalogue of Dental Products

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	<b>Page</b>
SAE Spark Erosion Machine	4/5
SAE Basic Equipment	6/7
SAE Erosion Accessories	8/9
Electrodes for Precious Metals / CoCrMo / Titanium	10/11
SAE Friction Pins / Laser Welding Wire	12/13
Telescopic Double Crowns with Controllable Friction	14/15
Dentures with Spark Erosion Attachments / Instructions of Use Customized Attachments	16/17
SAE T Attachment / SAE Special T Attachment	18/19
Chronological Process of a T Attachment Erosion	20/21
Assembly Instructions SAE Special T Attachment	22/23
SAE Double T Attachment	24/25
SAE Precision T 2000 Attachment	26-29
Instructions of Use SAE Swivel Latch Standard + Standard mini	30/31
SAE Swivel Latch Standard CoCrMo/Gold/Titanium	32/33
SAE Swivel Latch mini CoCrMo / Latch Gauge	34/35
Instructions of Use for SAE Swivel Latch Standard and Standard mini	35
SAE Module Swivel Latch / SAE Pin Attachment	36/37
Instructions of Use for SAE Module Swivel Latch	36/37
Assembly Instructions SAE Swivel Latch Standard + Standard mini	38/39
Assembly Instructions SAE Module Swivel Latch	40/41
Plastic Profiles / Tapered Cones / Diamond Polishing Pastes	42/43
SAE Model Centering System	44/45
Result of Comparative Study in QZ 03/2004 Implant Models with Constant Dimensional Stability	46/47
Precision Model Material / SAE Special Die Stone / Epoxy Resin and Hardener	48/49
SAE Precision Casting System for Crowns, Bridges and Single Cast / Disc with Fixing Pins	50/51
OKTA-Sil / Duplicating System	52/53
SAE Invest M2 / Metal Testing Cast / Octagonal Mould	54/55
Syringe / Dosage Dispenser / Wetting Agent / Bonding Agent / Activating Instruments	56/57
Alloys	58/59
Milling Burs	60/61
Certificates	62/63

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## → SAE Spark Erosion Machine

SAE Spark erosion machine EDM 2000

### Improved Performance: Spark Erosion Machine EDM 2000

- New digital generator technology
- Modernized processor controlled unit
- Technology data memory recall by product code
- Z axis unit
- Multi function display with error display
- Minimization of depth of surface roughness
- With PCI connector – laptop control possible
- Faster filling of tank
- Z axis “touch down zero point adjustment”
- Individual parameter adjustments possible

### Can be used for:

- SAE T attachment
- SAE double T attachment
- SAE swivel latch attachment
- SAE module swivel latch
- SAE pin attachment
- Friction pins
- Passivation of implant structures

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Basic Equipment for Spark Erosion Machine SAE EDM 2005/2000

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
## → Erosion Accessories

1. Chuck holder for securing chucks
- Chucks for securing copper electrodes, chuck capacity 0.5-6.0 mm**
2. Chuck 1.0-0.5 mm
3. Chuck 2.0-1.0 mm
4. Chuck 3.0-2.0 mm
5. Chuck 4.0-3.0 mm
6. Chuck 5.0-4.0 mm
7. Chuck 6.0-5.0 mm
8. Holder for securing suprastructure when eroding the implants
9. Universal electrode holder electrode holder for all T attachment erosions
10. Holder for centering system to achieve exact machine zero position
11. Holder for all latch attachment erosions
12. Depth gauge to adjust original depth for finishing erosion process
13. Electrode guide for erosion of latch axle
14. Angle bracket latch axle

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## → Accessories

- Dust cover (not shown)
- Filter for EDM 2000 and 2005 (not shown)
- Filter for SAE 100 and 1000 (not shown) from serial no. 1100 onwards
- Folded apron cover Width: 90 mm (not shown)
- Folded apron cover Width: 180 mm (not shown)
1. Cabinet

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Accessories for Spark Erosion Machine SAE EDM 2005 / 2000

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## → Model Mount with Angled Plate

1. Model mount for EDM 2005 and EDM 2000
2. Angled plate for model mount

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## → Attachment Device

1. Attachment Device for securing electrode band for making simultaneous slits in pins in the following sizes 1.6 + 1.7 + 1.8, slit size 0.11 mm

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## → Electrode Band

1. Electrode Band Copper  
Measurements:  
Strength 0.10 mm  
Width 10 mm  
Length 2000 mm  
for CoCrMo and gold

- Electrode Band Graphite (not shown)  
Measurements:  
Strength 0.10 mm  
Width 10 mm  
Length 45 mm  
for titanium

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## → Electrode Guide

1. Electrode guide for centering the round copper electrodes in the following sizes 0.6 - 1.7 mm Ø for all types of SAE spark erosion machines

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
Consumable Materials SAE Electrodes

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## → Copper Electrodes for Precious Metals/ CoCrMo /Titanium

Copper electrode Ø 0.6
Copper electrode Ø 0.7
Copper electrode Ø 0.8
Copper electrode Ø 0.85
Copper electrode Ø 0.9
Copper electrode Ø 1.0
Copper electrode Ø 1.1
Copper electrode Ø 1.2
Copper electrode Ø 1.4
Copper electrode Ø 1.5
Copper electrode Ø 1.6
Copper electrode Ø 1.7
Copper electrode Ø 1.9
Copper electrode Ø 2.2
Copper electrode Ø 2.7
Copper electrode Ø 2.9
Copper electrode Ø 4.0

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## → Carbide Electrodes for Titanium

Electrode Ø 0.80
Electrode Ø 1.50

### Tubular carbide electrodes for precious metals/CoCrMo / clover-leaf channell

Electrode Ø 0.80
Electrode Ø 0.90
Electrode Ø 1.50

#### Advantages of clover-leaf-channel tubular electrodes in tungsten carbide

- The specific clover-leaf form allows for improved whirl flushing. Additionally, this electrode is extremely robust and wears evenly, thus preventing distortion due to material burn off
- The channels are cylindrical throughout
- Ideal for blind holes since the floor of the hole stays flat
- Saves time

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
Consumable Materials SAE Attachments

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
## → Consumable Materials SAE T Attachment

1. SAE T attachment without screw
2. SAE T attachment with screw
- Replacement screw for SAE T attachment (not shown)
- | Composition (%) |      |     |      |      |       |       |  |
|-----------------|------|-----|------|------|-------|-------|--|
| Co              | Cr   | Mo  | Si   | Mn   | Fe    | N     |  |
| 64.2            | 28.0 | 5.9 | <1.0 | <1.0 | ≤0.75 | ≤0.25 |  |
3. Plastic pattern for attachment types:  
- SAE T attachment Standard without and with screw  
- SAE Special T attachment  
- SAE Precision T attachment
4. Copper electrode for approx. 15 erosions

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## → SAE Special T Attachment

- Special T attachment titanium incl. cone  
Height: 5.1 mm / 90°
1. Special T attachment titanium incl. cone  
Height: 3.0 mm / 90°
- Special T attachment titanium incl. cone  
Height: 5.1 mm / 45°
1. Special T attachment titanium incl. cone  
Height: 3.0 mm / 45°
2. Plastic pattern
3. Graphite electrode for SAE special T attachment, Length: 100 mm for 12 erosions
4. Copper electrode for SAE special T attachment for titanium, CoCrMo, precious metals, both sides can be used, Length: 100 mm for 25 erosions
- Replacement activating screw (not shown)
- Replacement screw for sealing cone (not shown)

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### Chronological Process of a T Attachment Erosion with Graphite and Copper Electrodes

1. Make a slit in the graphite electrode with a thin diamond disc approx. 5 – 6 mm below the erosion tip (see instruction no. 5). If working with copper, this step is not applicable.
2. Put electrode in holder with just approx. 20 mm protruding. Secure depth gauge to EDM 2000 / 2005 (order no. 30-2620).
3. Put model base with fixed model in originally chosen parallel position onto the work surface. Align surface of auxiliary part parallel to the shaped side of the electrode.
4. Now position the electrode using the XY adjustable axes\*\* over the erosion area so that there is buccally and palatinal to the matrix an even wall thickness of approx. 0.8 mm.
5. Remove primary part, lower electrode so that it almost touches the model. Move the electrode optically to the tooth abutment and then back again until there is a gap of 0.5 mm.
6. Determine height of auxiliary part (e.g. height or length is 5.3 mm, i.e. subsequent erosion depth is 4.8 mm), since remaining thickness must always be 0.5 mm.
7. Set up the primary part again and fix it using clamps (order no. 20-1018).
8. Position flushing hose at just under 45° so that it flushes away the burn-out from the erosion object into the tank.
9. Programme the generator to L 2 (gold LS 1): EDM 2000 / 2005 see technical data sheets. Start erosion process and allow to erode for a few seconds. Stop erosion process and make depth adjustment – as described in no. 6 (wear: graphite + 10 % / copper + 5 %). Start again and let it run.
10. Once process is complete, raise electrode. Carefully break off the part with the slit (if working with copper, exchange for electrode with sharp-edged front surface). Change the electrode once. Two erosion processes: rough, then smooth – use EDM 2000/2005 depth gauge.
11. Programme the generator to L 2 (gold LS 1) – as described in no. 9. Bring electrode optically to the same height as the auxiliary part and adjust the depth again which has already been identified\*. Briefly erode and then check the shape as to whether any small corrections must be made via the XY adjustment.  
\*Step not applicable with EDM 2000/2005 machine
12. After the erosion process, remove the primary part from the model, rinse and lightly blast with glass pearls. Adjust male component.

\*\* If an XY adjustment to the machine is not possible, then the adjustments have to be made manually.

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Consumable Materials SAE Double T Attachment

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


### → Double T Attachment

- 1. Double T attachment  
CoCrMo
- 2. Double T attachment  
Gold
- 3. Double T attachment  
Titanium

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### → Electrodes for Double T Attachment


- 1. Graphite electrode  
Double T attachment
- 2. Copper electrode  
Double T attachment

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
1. Complete kit 45° Titanium

2. Complete kit 45° Gold

Contents: 2 male parts complete assembled, 45°,  
2 plastic patterns, green

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## → Consumable Materials

1. Plastic pattern for Precision T 2000

2. Copper electrode for approx. 15 erosions

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1. Complete kit 45° Titanium

2. Complete kit 45° Gold

Contents: 2 male parts complete assembled, 45°,  
2 HSL matrices  
Retention sleeve in titanium

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**30-3935** 1. Matrize HSL

**30-3936** 2. Parallel-Fixator

## → Matrix HSL/ Paralleling Mandrel for Matrix

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1. Male part complete Titanium 45°

2. Male part complete Gold\* 45°

\* Female retention sleeve, screws and cone made of titanium

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**30-3910** 1. Patrizie Titan 45°

**30-3930** 2. Patrizie Gold 45°

## → Individual male parts Titanium/Gold

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### Swivel Latch Standard + Standard mini

#### Instructions of use

- Fig. 1: Position the 2° plastic profile using the parallelometer.
- Fig. 2: Finish off milling the cast tapered cone to 2°.
- Fig. 3: Position the prefabricated plastic latch onto the duplicate model to show necessary volume required for subsequent spark erosion.
- Fig. 4: Position the model to eroding position and finely adjust the primary structure at 90° angle to the occlusal or oral surface.
- Fig. 5: Selection of SAE electrodes for the necessary CSP and swivel latch attachment fit.  
Depth adjustment for electrode:  
The depth of 5.8 mm is to be regarded as a fixed dimension and must be set on the digital display at the beginning of the erosion process.
- Fig. 6: SAE friction pins and latch ready to be assembled.
- Fig. 7: Assembly of friction pins and swivel latch complete, specifically for unilateral dentures.

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## → SAE Model Centering System

- |   |
|---|
| 1. Machine holder for centering system                        |
| 2. Upper centering device                                     |
| 3. Centering pilot pin Parallelometer shaft thickness 2.34 mm |
| 4. Plastic centering mandrel                                  |

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### SAE Model Centering System

The model centering system was developed by SAE to determine the final 0 position (correct path of insertion) for the preliminary wax-up and for all subsequent steps when producing combined/fixed and removable prostheses. The model position, secured by the model mount, is fixed by centering the plastic mandrel with the help of the centering pilot pin. As such, the position cannot be changed during the following procedures.

#### Areas of use:

- Model preparation
- Milling and transfer work
- Erosion work

#### Working Procedures

1. Align working model or master model with centering pilot pin (order no. 60-1001) in the parallelometer to determine the 0 position.
2. Then mill a blind hole into the model to accommodate plastic centering mandrel (order no. 60-1003) in 0 position. This is achieved by using a vertical milling machine which is equipped with a star drill bit.
3. Positioning the plastic centering mandrel in the 0 position is done in the parallelometer.
4. The plastic centering mandrel is pushed onto the centering pin of the parallelometer and then inserted into the milled blind hole using a vertical movement.
5. The mandrel is secured in the model using pattern resin. The precision fit of the plastic centering mandrel guarantees the transfer of the 0 position for the working procedures – modelling, milling and erosion machining.
6. Align with already determined 0 position:

By lowering the centering pilot pin to fit into the plastic centering mandrel, the mounted working or master model automatically attains the 0 position.

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## → Duplicating System

1. Duplicating base plate  
plus plasticine  
– reusable –

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**70-1004** 2. Dublierständer

## → Stabilizing Plate and Duplicating Stand

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


## → SAE Invest M2

- SAE Invest M2  
Package / Contents:
1. 20 x 200 g bag for model
  2. 20 x 500 g bag for outer mould and 1.8 litres special liquid

SAE Invest M2 is a phosphate bonded investment for model casts. SAE Invest M2 is mixed with a special liquid which enables adjustable expansion. Secondary parts can therefore be cast as well with an accurate fit using the single cast method.

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## → Metal Testing Cast

1. Metal testing  
cast for adjusting  
expansion of the  
SAE Invest M2

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- 70-4003** per Packung = 100 Stück

## → Octagonal Mould

1. Octagonal mould  
Package includes:  
1 side segment with funnel  
access,  
3 side segments,  
1 base plate,  
5 funnels

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## → Syringe and Dosage Dispenser

1. Syringe, 60 ml  
 2. Duplicating dosage dispenser with pump system for silicone materials Capacity 3.5 l

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## → Wetting Agent, Bonding Agent and Diluting Agent

1. Wetting Agent Eliminates surface tension of model before filling silicone into duplicating mould, 500 ml  
 2. Bonding Agent Increases wax bonding of Invest M2 when using prefabricated wax patterns, 30 ml  
 3. Bonding Agent, as 70-1012, 100 ml  
 4. Diluting Agent, 100 ml

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## → Special Alloy OKTA-C

OKTA-C – nickel free  
Special alloy  
for crowns and bridges,  
for ceramic veneer

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## → Special Alloy OKTA-M VS

OKTA-M VS – nickel free  
Special alloy  
for claspless casts  
for precision  
single cast method  
for combined restorations

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## → Screw Connections

- 1. Cylinder head screw titanium for threaded sleeve, small, M 1.6 mm
- 2. Grub screw titanium for threaded sleeve, small, M 1.6 mm
- 3. Threaded sleeve titanium, small, M 1.6 mm
- 4. Threaded sleeve CoCrMo, small, M 1.6 mm
- 5. Cylinder head screw large, titanium, M 2.0 mm
- 6. Grub screw titanium for titanium threaded sleeve, large, M 2.0 mm
- 7. Grub screw gold for gold threaded sleeve, large, M 2.0 mm
- 8. Threaded sleeve CrCoMo, large, M 2.0 mm
- 9. Threaded sleeve titanium, large, M 2.0 mm
- 10. Threaded sleeve gold, large, M 2.0 mm
- 11. Driver
- 12. Driver - hand driver

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**\*These system parts are used only for making dental prostheses!**

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

## → Handpiece Milling Burs

- 1. Handpiece milling bur, high speed, carbide, pointed, oval, small
- 2. Handpiece milling bur, high speed, carbide, pointed, oval, large
- 3. Handpiece milling bur, high speed, carbide, pointed, small
- 4. Handpiece milling bur, high speed, carbide, pointed, large
- 5. High performance handpiece milling bur, carbide, conical shape

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## → ZR Grinders

Grind zirconium oxide easily with:

- 1. ZR grinder for friction pin Ø 1.0 rough (no colour)
- 2. ZR grinder for friction pin Ø 1.0 medium (red)
- 3. ZR grinder for friction pin Ø 1.0 fine (yellow)

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## → 2° Milling Burs

- 1. End mill 2.35 mm for screw channel
- 2. Screw based milling bur

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## → Guide bushing and lock nut for passivation of screw channel

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## Certificate

**mdc medical device certification GmbH**  
 certifies that

**SAE DENTAL VERTRIEBS GMBH  
 - INTERNATIONAL -**  
**Langener Landstraße 173  
 27580 Bremerhaven  
 Germany**  
 for the scope

development, procuring, manufacturing, final inspection and sales of  
 dental precision fastenings, dental alloys, implant accessories and  
 spark erosion machines and accessories for the dental field

has introduced and applies a

### Quality Management System

The mdc audit has proven that this quality management system  
 meets all requirements of the following standard

**EN ISO 13485**

Medical devices – Quality management systems –  
 Requirements for regulatory purposes

EN ISO 13485:2012 • AC 2012 • ISO 13485:2003 • Csk 1:2009

Valid from	2014-03-27
Valid until	2019-03-26
Registration no.	D1083300007
Report no.	E-D10833_2014-03-19
Stuttgart	2014-03-19

*J. Borch*  
 Head of Certification Body

mdc medical device certification GmbH  
 Kriegerstraße 6  
 D-70191 Stuttgart, Germany  
 Phone: +49 (0)711-253597-0  
 Fax: +49 (0)711-253597-10  
 Internet: <http://www.mdc-cd.de>

## EC Certificate

**mdc medical device certification GmbH**  
 Notified Body 0483  
 herewith certifies that

**SAE DENTAL VERTRIEBS GMBH  
 - INTERNATIONAL -**  
**Langener Landstraße 173  
 27580 Bremerhaven  
 Germany**  
 for the scope

dental precision fastenings and dental alloys

has introduced and applies a

### Quality System

for the manufacture of the products concerned and carries out a  
 final inspection as specified in Annex V, Section 3.

The mdc audit has proven that this quality system  
 meets all requirements according to

**Annex V – Section 3  
 of the Council Directive 93/42/EEC**

of 14 June 1993 concerning medical devices.  
 The surveillance will be held as specified in Annex V, Section 4.

Valid from	2014-03-27
Valid until	2019-03-26
Registration no.	D1083300008
Report no.	E-D10833_2014-03-19
Stuttgart	2014-03-19

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