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**IMPLANT-retained suprastructures
passivated to a stress free fit through spark erosion**

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

See: Quintessenz Journal of Dental Technology
Volume 2 / Number 2 / April 2004, page 158-181
Edentulous and partially dentate patients:
Quality assurance through interdisciplinary
cooperation

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See in Korean: QDT 5 / Volume 28 / 2003

See in Polish: Quintessence Periodontologia Implanty
2 / May 2004

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See: EDI-Journal 02/2006
Passivation of implant superstructures using the
spark-erosion technique

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

**One outer screw is fastened
before spark erosion passivation of the structure
Misfit!**

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**One center screw is fastened
before spark erosion passivation of the structure
Misfit!**

SAE-SECOTEC / Sheffield-Test

The spark erosion method of correcting inherent inaccuracies of conventional prosthetic techniques when manufacturing implant-retained superstructures.

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*All screws are fastened
before spark erosion passivation of the structure
Stress transfer to the implants*

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*One outer screw is fastened
after spark erosion passivation of the structure
stress-free fit of the structure on the implants*

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

The SAE Spark Erosion Method for Achieving a Passive Fit of Implant Framework

The conventional methods of casting framework in the dental laboratory will never lead to a passive fit of the sub and suprastructures situated on osseointegrated implants.




However, to comply with the demands of dentistry to achieve a passive fit of the implant related framework, many variables are considered. The cast suprastructure is segmented and reconnected, either by soldering or welding. This procedure leads to new inaccuracies, if only to an insignificant extent. But when the ceramic veneering is applied to the multispanded restoration, the framework shows a stress producing potential due to the shrinkage of the ceramic during firing. Again, we find an inaccurate fitting fixture.

These inaccuracies resulting from conventional prosthetic techniques can be corrected after completion of the suprastructure by applying the SAE Secotec Spark Erosion Method.

With the SAE Secotec Spark Erosion Technique, short circuit impulses in an electrical field burst out particles from the metallic surface of the suprastructure, starting from the areas that first come into contact with the lab implants and continuing until all contours of the framework are equally spark eroded, thus achieving uniform contact.

A special model is needed to provide electrical contact during the spark erosion procedure. The model shells furnish the opportunity to exchange the implant replicas for identical electrodes assuring an even flow of electricity and a sufficient production of sparks.

The Secotec System is applicable to all conductible metals and alloys and to acrylic or ceramic veneered suprastructures.

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What is Spark Erosion

The Spark Erosion is a method for harming design of electric-conductive metals by controlled electric electricity impulses between tool electrode (implant electrode) and working piece (implant meso or suprastructure) in presence of a non-conducting (dielectrical) liquid.

The erosion is compared with the machining procedure (cutting – drilling) an additional procedure for working metals, but however with the essential difference that the rotation of the tool (drill cutter) is dropped, and there is created the opportunity by the fixed electrode in the vertical propulsion especially in case of restorations on implants to obtain utmost fits with individual design.


For 20 years the Spark Erosion System has been used in the Dental Technology in the area of the combined prosthodontia. In the implant prosthodontia the procedure has been used for several years for passivation of the fit with regard to meso and suprastructures on implants.

The contraction caused by the cast of implant meso and suprastructures can quickly be removed by the spark erosion system of the contact areas causing tension, so that a passive sit of the meso and suprastructures on the implants can be obtained.

The principle of the Spark Erosion


Working piece (meso or suprastructure) and electrode are turned in the spark erosion machine into working position (illustration 1 and 2), so that they do not touch each other. It rests a gap which is filled with an isolating liquid (dielectricum). The working on it is effected in a container. Working piece and electrode are connected by a current wire with a direct-current source. The electric cables are interrupted by a switch. If it is closed, there are caused, by short-circuit, impulses between working piece and electrode, which emit sparking somersaults from the electrode to the working piece. These sparking somersaults cause the harming of the metal on the working piece (structure) on the parties which hit first. The harming continues until all parties of the structure have full contact with the electrodes. The working piece sits fully circulative on the implant electrodes whose surface design corresponds exactly with the surface design of the implants or of the abutments.

The Secotec System under Use of the Spark Erosion

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
An electric field is developed between electrode and work piece by turning on the tension.

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
The electrically conducted particles concentrate in the field between electrode and work piece and form a bridge.

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
Suddenly a channel of discharge develops itself over this bridge.

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
With beginning discharge, the number of the electrically conducted particles increases like an avalanche, the electricity begins to flow.

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
A counter-pressure arises because of the viscosity of the dielectric which ties up the channel of discharge between electrode and working piece. The concentrated heat brings the surface areas to melting.

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
The temperature increases, the melting continues into the surface areas. High temperatures arise by the vaporization of the melt in short-term.

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The flow of electricity is interrupted by switching off the tension, the channel of discharge fails. This leads to sudden centrifugation and vaporization of the melt.

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There is spread fine electrode material in the dielectric at the end of the discharge. The surface areas of the working piece and the electrode have an ablation crater each.

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Purpose of the Dielectricum

Insulation:

The dielectricum has the important function to insulate the working piece from the electrode. The disruptive discharge must be effected at the smallest distance. Thereby the ablation capacity increases and the precision of reproduction.

Cooling:

Since the electric spark hits on the working piece with about 10,000°C, the dielectricum must cool electrode and working piece. An overheating of the electrode must be avoided, so that there does not arise an augmented high loss. The metal steams developed by the Spark Erosion System must be able to condense in the liquid, additionally the resulted ablation particles must be washed away from the erosion zone to avoid faults with regard to the process.

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
Figs. 1 + 1a

The individual tray holding the impression in Impregum with the impression copings.

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Figs. 2 + 2a

The lab implant replicas of a matching implant system are attached to the Secotec model shells and then screwed to the impression copings embedded in the impression. All tightening procedures are performed with a torque wrench (item no. 82-0521) and a counter wrench (item no. 82-0531), hereby complying with the specific torque – Ncm. – Copper wire (item no. 82-0500) is attached to each model shell, thus connecting them to form an electrical circuit.

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Fig. 3

An elastic die acrylic is poured. The implant replicas are completely covered with the soft acrylic and only the model shells are visible.

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Fig. 5

The entire area is slightly dusted with silver powder to insulate against the Epoxy Resin.

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Fig. 6

A partial pour is made with low contracting – 0.03 mm – die Epoxy Resin (item no. 40-1060).

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Fig. 7


After all of the excess wax is removed the remaining pour is completed in a type IV die stone.

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Model Structure

- 10 die model
- 10a removable elastic portion
- 12 admitting area for implant replica and electrode
- 12a model shell
- 14 threading path of model shell
- 16 copper wire contact
- 18 copper wire for conductivity (anode)

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



Fig. 8

- a) torque wrench
- b) driver
- c) implant replica
- d) model shell
- e) counter wrench

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
SAE Spark Erosion Machine

10-1105 SAE Spark erosion machine EDM 2005 **new**

Improved Performances: Spark Erosion Machine EDM 2005

- New digital generator technology
- Modernized processor controlled unit
- Technology data memory Call up by product code
- Z axis unit
- Multi function display with error display
- Minimization of the surface rough depth
- With PCI contact – control by laptop possible
- Faster filling of tank
- Z axis “touch down zero point adjustment”
- Individual parameter adjustments possible

Usable for:

- SAE T attachment
- SAE double T attachment
- SAE interlock attachment
- SAE swivel-latch attachment
- SAE swivel-latch attachment '95
- SAE module latch
- SAE pin attachment
- Friction pins
- Passivation of implant frameworks

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Erosion Accessories

1. Chuck holder to receive chucks
Chucks to receive copper electrodes, chuck capacity 0.5-6.0 mm
2. Chuck holder 1.0-0.5 mm
3. Chuck holder 2.0-1.0 mm
4. Chuck holder 3.0-2.0 mm
5. Chuck holder 4.0-3.0 mm
6. Chuck holder 5.0-4.0 mm
7. Chuck holder 6.0-5.0 mm
8. Holder to engage superstructure for erosive machining with implant abutments
9. Universal electrode holder electrode holder for all T attachment erosions
10. Holder for centering system to achieve exact 0 position in machine
11. Holder for all latch attachment erosions
12. Depth gauge to enable the adjustment of original depth for finishing erosive process
13. Electrode guide for latch axis
14. Angle holder latch axis

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Accessories

15. Model mount for EDM 2005 and EDM 2000
16. Angled plate for model mount
17. Securing clamps for model (1 pair)
18. Clamp
Rubber circulars 0.5 - 2.5 mm for EDM 2005 and EDM 2000 (not shown)
Dielectric fluid (not shown)
Sleeves for electrode guide for 30-1602 (not shown)

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Ankylos

1. model shell standard
2. implant replica
3. implant electrode 4,45 Ø
4. plastic burn-out cylinder
5. screw ❶
6. waxing screw
7. impression screw
8. imprint cap
9. driver implant replica
10. driver electrode and screw
11. driver waxing screw

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Astra 20°

1. model shell standard
2. implant replica
3. implant electrode
4. plastic burn-out cylinder
5. screw ❶
6. driver implant replica
7. driver electrode and screw

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Bego-Semados

1. model shell standard
2. implant replica
3. implant electrode
4. plastic burn-out cylinder
5. screw ❶
6. waxing screw
7. driver implant replica
8. driver electrode and screw

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


3i Certain Ø 4.0 mm

1. model shell standard
2. implant replica Ø 4.0 mm
3. implant electrode Ø 4.0 mm
4. plastic burn-out cylinder Ø 4.0 mm
5. screw hexagon ①
6. screw square ①
7. driver implant replica
8. driver electrode and screw
9. driver screw square

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3i Certain Ø 5.0 mm

1. model shell standard
2. implant replica Ø 5.0 mm
3. implant electrode Ø 5.0 mm
4. plastic burn-out cylinder Ø 5.0 mm
5. screw hexagon ①
6. screw square ①
7. driver implant replica
8. driver electrode and screw
9. driver screw square

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3i Certain Ø 6.0 mm

1. model shell standard
2. implant replica Ø 6.0 mm
3. implant electrode Ø 6.0 mm
4. plastic burn-out cylinder Ø 6.0 mm
5. screw hexagon ①
6. screw square ①
7. driver implant replica
8. driver electrode and screw
9. driver screw square

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**3i Conical Abutment
WCA 51-55**

1. model shell standard
2. implant replica
3. implant electrode
4. plastic burn-out cylinder
5. screw ①
6. driver implant replica
7. driver electrode and screw

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▶ Nobel Biocare Replace Select Ø 6.0 mm

1. model shell standard
2. implant replica
3. implant electrode
4. plastic burn-out cylinder
5. screw 2,0 ❶
6. driver implant replica
7. driver electrode and screw

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This system phases out and is replaced with Nobel Biocare Multi Unit. Identical with Mirus-Cone Ø 4.8 RP.

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▶ Nobel Biocare UCLA

1. model shell
2. implant replica
3. implant electrode
4. plastic burn-out cylinder
5. screw ❶
6. direction screw
7. driver electrode and screw

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N15 ▶ SAE-Secotec für UCLA-WP (Wide Platform)

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Schütz-Impla

1. model shell standard
2. implant replica
3. implant electrode 4,45 Ø
4. plastic burn-out cylinder
5. screw ❶
6. abutment 3.3 Ø/hight 1.5+3.0*
7. abutment 4.2 Ø/hight 1.5+3.0*
8. impression screw
9. imprint cap
10. driver implant replica
11. driver electrode and screw
12. driver abutment

* original system component/Abutment

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S10 ▶ SAE-Secotec für **straumann SynOcta WN** (Wide Neck)

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Straumann SynOcta

1. model shell standard
2. implant replica
- * SynOcta secondary part
(original system component/Abutment)
3. implant electrode
4. plastic burn-out cylinder
5. screw ❶
6. driver implant replica
7. driver SCS screw
8. driver electrode and screw

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Straumann SynOcta 1.5

1. model shell standard
2. implant replica
- * SynOcta 1.5 secondary part
(original system component/Abutment)
3. implant electrode
4. plastic burn-out cylinder
5. screw ❶
6. driver implant replica
7. driver SCS screw
8. driver electrode and screw

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Straumann SynOcta

1. model shell standard
2. implant replica
- * SynOcta secondary part
(original system component/Abutment)
3. implant electrode
4. plastic burn-out cylinder
5. driver implant replica
6. driver SCS screw
7. driver electrode and screw

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Straumann SynOcta WN

1. model shell standard
2. implant replica
- * SynOcta WN secondary part
(original system component/Abutment)
3. implant electrode
4. plastic burn-out cylinder
5. screw ❶
6. driver implant replica
7. driver SCS screw
8. driver electrode and screw

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▶ Straumann Narrow-Neck-Implantat

1. model shell standard
2. implant replica
3. implant electrode
4. plastic burn-out cylinder
5. screw ①
6. driver
7. driver electrode and screw

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▶ SPI Thommen Ø 3.5

1. model shell standard
2. implant replica
3. implant electrode
4. plastic burn-out cylinder
5. screw M 1.4 ①
6. driver implant replica
7. driver electrode and screw

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▶ SPI Thommen Ø 4,0

1. model shell standard
2. implant replica
3. implant electrode
4. plastic burn-out cylinder
5. screw M 1.6 ①
6. driver implant replica
7. driver electrode and screw

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

* Originalsystemteil/Abutment

SPI Thommen Ø 5.0

1. model shell standard
2. implant replica
3. implant electrode
4. plastic burn-out cylinder
5. screw M 1.6 ❶
6. driver implant replica
7. driver electrode and screw

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Zimmer Spline: Ø 3.25

1. model shell
2. implant electrode GH 1
3. implant electrode GH 2
4. implant electrode GH 3
5. plastic burn-out cylinder
6. screw ❶
7. driver electrode and screw

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Transparent model

1. transparent model holding
4 implant replicas
2. transparent model holding
4 implant replicas and dolder bar
in titanium – laser welded
3. Complete transparent model with dentures

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Implant replica

1. implant replica
2. conical cap
3. screw
transparent model compound per kg –
not pictured

**Implant abutment replicas are intended
only for demonstration models and not
for clinical cases.**

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Titanium cylinder

1. titanium cylinder inner Ø 3.1
for pilot drill
2. titanium cylinder
inner Ø 2.00 for Lindemann miller
length 8 mm, fits into 82-0497
3. titanium cylinder inner Ø 2,4

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Mixing Ratio 1 to 6:

volume: hardener 1 ml: Resin 6 ml
weight: hardener 1 g: Resin 10 g
(see also QZ 3/2004, pages 238 - 250)

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Information on Implant-retained Restorations

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Information on Implant-retained Restorations – Concerning Screws ①

Implant-related screws of the Secotec System manufactured by SAE are for laboratory use only. To clinically secure implant-retained mesostructures and superstructures please use the original systems screws – no SAE Secotec screws – risk of breakage!

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