


SAE Horizontal Centrifuge for the consistently shaped master model

Made in Germany



Sectioned cast



Master model




Implant model

SAE Horizontal Centrifuge for the bubble-free epoxy resin model

Technical data:

- Measurements:
W: 62 cm D: 51 cm H: 35 cm (closed) / 55 cm (open)
- Weight: 45 kg
- Powerful DC motor with 240 W
- Revolution centrifuge adjustable from 0 – 500 rpm
- Can be used worldwide thanks to the efficient (up to 87%) switch-mode power supply
- Input voltage:
85 – 264 VAC @ 47 – 63 HZT or 120 – 370 VDC
- Maximum input power: 336 W

SAE Horizontal Centrifuge: Art. No. 55-0780



SAE DENTAL VERTRIEBS GMBH

Langener Landstraße 173 · D-27580 Bremerhaven
Tel.: +49 (0)471 9 84 87 45 · Fax: +49 (0)471 9 84 87 44
E-Mail: info@sae-dental.de · Internet: www.sae-dental.de